



# Purus Gas Systems

Top: Danielle Hicks, Cost Accountant

Bottom: David Smith, Chief Scientist, Restek Performance Coatings



for **more** info**Questions about which carrier gas purifier to use?**

Call 800-356-1688 or 814-353-1300, ext. 4, or contact your Restek representative to discuss your application with our technical service chemists.

**Why do I need to use traps and where should I install them?**

Carrier gas must contain less than 1ppm of oxygen, water vapor, or any other trace contaminant, to prevent column degradation, shortened column lifetime, and increased stationary phase bleed. Contaminants cause ghost peaks to appear during temperature programming and degrade the validity of analytical data. Make-up gas also should be contaminant-free, or baseline fluctuations and excessive detector noise can occur; detector gases should be free of water and hydrocarbons, or excessive baseline noise can result. Gas purifiers remove these contaminants from gas sources, thereby improving system performance.

**Moisture Removal**

Moisture in carrier gas lines will prematurely degrade oxygen and hydrocarbon traps and increase detector noise (particularly with ECDs). As a precaution, we highly recommend installing a moisture trap before the hydrocarbon and oxygen traps on all carrier gas lines. Our favorite trap is the Super-Clean™ Ultra-High Capacity Moisture Filter (cat.# 22028, pg. 220).

**Hydrocarbon Removal**

Use a hydrocarbon trap if your gas has a potential source of hydrocarbon contaminants (e.g., an oil pump in an air compressor) or if you suspect you are observing carrier gas ghost peaks. Install the hydrocarbon trap after the moisture trap, to prevent moisture from degrading the hydrocarbon-trapping ability of the activated carbon in the hydrocarbon trap. We recommend the Super-Clean™ Ultra-High Capacity Hydrocarbon Filter (cat.# 22030, pg. 220).

**Oxygen Removal**

Oxygen is a column killer. It is present even in UHP gases, as minute leaks at fittings allow oxygen to influx against the concentration gradient. There are many choices for oxygen removal—the Super-Clean™ Ultra-High Capacity Oxygen Filter (cat.# 22029, pg. 220) is popular with Restek chemists. Because oxygen can enter a gas line at any fitting, the oxygen trap should be the last connection before the gas line enters the chromatograph.

**tech tip****Carrier Gas Purity**

Carrier gas should contain less than 1ppm of oxygen, moisture, or other trace contaminants, to prevent column degradation and to decrease stationary phase bleed.

The cost of high-purity gases and carrier gas line purifiers is offset by longer column lifetimes, more accurate data, and less GC maintenance.

did you **know?****Trap replacement made simple!**

Try the Super-Clean™ Triple Filter Carrier Gas Cleaning Kit (cat.# 22019, pg. 220)—it removes moisture, hydrocarbons, and oxygen in one easy-to-change, economical cartridge.

# Elements of Gas System Design



## Purus™ Gas Systems

We know how important high-purity gas is to the success of your analysis. To provide you with the products and services you need to achieve a first-quality gas delivery system, Restek has created Purus™ Gas Systems. Purus™ Gas Systems are application-specific assemblies of the finest products and technical expertise for a systematic approach to gas handling. "Systematic" means that every product needed to deliver research purity gas is available from one source, at the level of quality you expect from Restek. We will not only supply products, we will work with you to design the best gas system for your application.

GC ACCESSORIES | PURUS GAS SYSTEMS



### Gas Source

### Source Pressure Control

### Tubing, Fittings, & Tools


### Preliminary Gas Purification

#### Low-Pressure Gases





#### High-Pressure Gases

##### Direct Pressure Control



##### Point-of-Use Pressure Control



We provide the following gas system products:


- Restek gas purifiers provide cost-effective gas purity assurance solutions.
- Restek stainless steel and copper tubing, cleaned and ready to use.
- Parker A-Lok® tube fittings consistently deliver high-quality performance.
- Extensive line of hand tools to make your work easier and faster.

Restek's Technical Service Team, 800-356-1688 or 814-353-1300, ext. 4, or your Restek representative, can answer your questions and provide technical support literature or system-design advice. From the gas source to your point of use, we offer the products and services that ensure the purity of your gas.

## Line Pressure Control

## Final Gas Purification

## GC



# Thermal Gas Purifier, Restek Gas Management System

## Thermal Gas Purifier\*

- Removes oxygen, water, carbon monoxide, carbon dioxide, hydrocarbons (except methane) to ppb levels—pure enough for MS.
- Packed with reactor-grade, pure getter material for maximum efficiency and no contamination.
- Welded end-fittings on getter tubes eliminate leaks.
- Each tube has 12L oxygen and 35L water vapor capacity.
- Maximum flow: 1 liter/minute.




Single-Tube Model



Dual-Tube Model

The getter material in Restek's re-engineered line of thermal gas purifiers reacts chemically with impurities in the carrier gas stream. Because the reaction is non-reversible, there is no possibility of contaminants breaking through the thermal gas purifier.

Gas purification is very economical when you use a thermal gas purifier. Getter tubes normally require changing only once per year; heavy use or very impure feed gas may require more frequent getter tube replacement.



welded end fittings for leak-free connection

Dimensions: 13" x 1/2"  
(33 x 1.3 cm)

Restek Single-Tube Thermal Gas Purifier, 110 Volt (1 tube included)**	qty.	cat.#
1/8" Fittings	ea.	21496
1/4" Fittings	ea.	21497
Restek Dual-Tube Thermal Gas Purifier, 110 Volt (2 tubes included)**	qty.	cat.#
1/8" Fittings	ea.	21498
1/4" Fittings	ea.	21499
Replacement Getter Tube	qty.	cat.#
1/8" Fittings (Similar to Supelco part# 2-2396)	ea.	21661
1/4" Fittings (Similar to Supelco part# 2-2398)	ea.	21660

\*Use with helium or nitrogen carrier gas only. Do not use with hydrogen, oxygen, or air—a violent reaction and/or fire will occur.

\*\*Single-tube model dimensions: 15" x 7" x 6" (38 x 17.8 x 15.2 cm). Dual-tube model dimensions: 15" x 10" x 6" (38 x 25.4 x 15.2 cm).



## Pressure Gauge Kit

- Use an in-line pressure gauge to indicate when the Thermal Gas Purifier getter tube should be replaced.
- Includes 1/8" tee and 0–100 psi gauge.

Always know the status of the Thermal Gas Purifier getter tube: change tube when pressure drops below 10psi.

Description	qty.	cat.#
In-line Pressure Gauge Kit for Thermal Gas Purifiers	kit	21657


## Restek Gas Management System

- Removes moisture, hydrocarbons, and oxygen from carrier gas, extending column lifetime.
- Produces high-purity carrier gas for most applications.
- Includes one each: moisture, hydrocarbon, and indicating oxygen trap.
- Replacing traps is safe and easy.
- Maximum flow: 1 liter/minute.



Dimensions: 12" x 14" x 3"  
(30.5 x 35.6 x 7.6 cm)

trap configuration →



## did you know?

The Restek Gas Management System removes water vapor (to 10ppb), hydrocarbons (to 0.1ppm), and oxygen (to less than 0.1ppm) with three traps housed in one unit.

Description	qty.	cat.#
Restek Gas Management System (includes fittings for 1/8" and 1/4" gas line)	ea.	21999
Replacement Traps	qty.	cat.#
High-Capacity Moisture Trap, 1/8" Fittings	ea.	21997
Capillary-Grade Hydrocarbon Trap, 1/8" Fittings	ea.	21991
Indicating Oxygen Trap, 1/8" Fittings	ea.	22010

**Super-Clean™ Gas Filters**

- High-purity output ensures 99.9999% pure gas.
- “Quick connect” fittings for easy, leak-tight cartridge changes.
- Glass inside to prevent diffusion; plastic outside for safety.

**Designed for fast, simple cartridge changing**

Cartridge systems make changing gas filters quick and easy, and Super-Clean™ gas filters are the latest in cartridge-style gas filtration. A baseplate allows cartridges to be exchanged without introducing oxygen. Spring-loaded check valves seal when a filter is removed and open only when a new filter has been locked in place. There is no longer a need for loosening and tightening fittings every time a trap is changed, and your system will not become contaminated during the process.

With available 2- or 3-position baseplates, you can purify all GC gas streams at one location. Figure 1 shows some possible filter cartridge combinations using these baseplates. Any combination is possible because any Super-Clean™ filter cartridge can be used with any baseplate.

**High-purity output for improved sensitivity (Table I)**

The Triple Filter cartridge (cat.# 22020) is ideal for purifying carrier gas. This trap contains oxygen, moisture, and hydrocarbon scrubbers in one cartridge. The purity of your carrier gas after flowing through the Triple Filter is better than six-9s (99.9999% pure), which is ideal for sensitive mass spectrometry (MS) or ECD analyses, and for protecting your analytical columns against damage from contaminated carrier gas.


The Fuel Gas Filter cartridge (cat.# 22022) is perfect for purifying flame ionization detector (FID) fuel gases, removing both moisture and hydrocarbons. Using the Fuel Gas Filter for FID hydrogen and air will produce a stable baseline, improving overall reproducibility and sensitivity.



**did you know?**

All Super-Clean™ filter cartridges (except hydrocarbon filter cat.# 22030) feature easy-to-read indicators. The indicator code is shown on every trap so there is no confusion about when to replace it.

**Figure 1** Filter cartridges can be configured for different applications.



**Table I** Each Super-Clean™ filter provides high-purity outlet gas.

Type of Filter	Outlet Gas Quality (%)	Maximum Pressure/ Maximum Flow Rates	Use for:	Indicator Color Change	Capacity			Estimated Lifetime (years)
					H <sub>2</sub> O (g)	O <sub>2</sub> (mL)	Hydrocarbons (g)	
Moisture cat.# 22028	>99.9999	11 bar 159psi/ 7 L/min.	Inert carrier gas Air Hydrogen	Yellow to clear	7.2	—	—	>2
Oxygen cat.# 22029	>99.9999	11 bar 159psi/ 7 L/min.	Inert carrier gas	Green to grey	NA	1000	—	>2
Hydrocarbons cat.# 22030	>99.9999	11 bar 159psi/ 7L/min.	Inert carrier gas Air Hydrogen	No indicator	NA	—	12 <sup>3</sup>	>2
Fuel Gas <sup>1</sup> cat.# 22022	>99.9999	11 bar 159psi/ 7L/min.	Inert carrier gas Air Hydrogen	Yellow to clear	3.5	—	24 <sup>3</sup>	>1.5
Triple <sup>2</sup> cat.# 22020	>99.9999	11 bar 159psi/ 7L/min.	Inert carrier gas	Yellow to clear Green to grey	1.8	500	4 <sup>3</sup>	>1
Helium Specific <sup>2</sup> cat.# 21982	>99.9999	11 bar 159psi/ 7L/min.	Helium	Yellow to clear Green to grey	1.8	500	—	>1

<sup>1</sup>Removes hydrocarbons, moisture.

<sup>2</sup>Removes hydrocarbons, moisture, oxygen.

<sup>3</sup>As a butane


**tech tip**

**Oxygen and Moisture Traps**

We highly recommend oxygen and moisture traps for make-up gas when operating sensitive detectors such as electron capture detectors (ECD). The hydrogen reaction gas for sensitive electrolytic conductivity detectors (ELCD) also requires a hydrocarbon trap, to remove trace impurities.

**Super-Clean™ filters are listed on pages 220 and 221.**

## Super-Clean™ Filters



All traps measure: 10<sup>5</sup>/<sub>8</sub>" x 1<sup>3</sup>/<sub>4</sub>"  
(27 x 4.4 cm)

Each baseplate unit measures:  
4" x 4" x 1<sup>7</sup>/<sub>8</sub>"  
(10.2 x 10.2 x 4.8 cm)

### Super-Clean™ Filter and Baseplate Kits

Description	qty.	cat.#
Carrier Gas Cleaning Kit (includes mounting baseplate, 1/8" inlet/outlet fittings, and oxygen/moisture/hydrocarbon Triple Filter)	kit	22019
Fuel Gas Purification Kit (includes mounting baseplate, 1/8" inlet/outlet fittings, and hydrocarbon/moisture Fuel Gas Filter)	kit	22021


### Replacement Filters



Description	qty.	cat.#
Replacement Triple Filter (removes oxygen, moisture and hydrocarbons)	ea.	22020
Replacement Fuel Gas Filter (removes moisture and hydrocarbons)	ea.	22022

### Filter Bundle Kit

Kit includes two Fuel Gas Filters for FID fuel gases and one Triple Filter for carrier gas.



Description	qty.	cat.#
Filter Bundle Kit	kit	22031

### Super-Clean™ Ultra-High Capacity Filters



Description	qty.	cat.#
Ultra-High Capacity Hydrocarbon Filter	ea.	22030
Ultra-High Capacity Moisture Filter	ea.	22028
Ultra-High Capacity Oxygen Filter	ea.	22029

### Helium-Specific Super-Clean™ Filter and Kit

- Specifically designed for purification of helium in GC/MS systems.
- Traps are packed and conditioned using helium.
- Uses standard single-position baseplate.



Description	qty.	cat.#
Helium-Specific Carrier Gas Cleaning Kit (includes mounting baseplate, 1/8" inlet/outlet fittings, and oxygen/moisture/hydrocarbon Helium-Specific Filter)	kit	21983
Replacement Helium-Specific Filter (removes oxygen, moisture and hydrocarbons)	ea.	21982

### Baseplates

Standard baseplate fittings are 1/8". To adapt to 1/4", order 1/8" to 1/4" tube-end unions (cat. # 21833, below).



Description	qty.	cat.#
Single-Position Baseplate	ea.	22025
2-Position Baseplate	ea.	22026
3-Position Baseplate	ea.	22027

### Wall Mounting Bracket

Baseplates can be mounted by using screws and the mounting holes on the baseplate, or by using this optional wall mounting bracket.



Description	qty.	cat.#
Wall Mounting Bracket for Super-Clean™ Baseplates	ea.	21984

### Replacement O-Rings for Cartridge Baseplates

Pack includes 10 large O-rings and 10 small O-rings.



Description	qty.	cat.#
Replacement O-Rings for Cartridge Baseplates	20-pk.	22023

### 1/8-Inch to 1/4-Inch Tube-End Unions

To adapt 1/8" Super-Clean™ baseplate fittings to 1/4", use 1/8" to 1/4" tube-end unions.



Description	qty.	cat.#
1/8" to 1/4" Tube-End Unions	5-pk.	21833

## Super-Clean™ Gas Trapping System for LC/MS

### Quick-change cartridge system for removing hydrocarbon impurities from nitrogen

- Changing filters is quick and easy.
- Up to 20L of hydrocarbon-free nitrogen per minute.

The Super-Clean™ Gas Trapping System is the latest technology in cartridge-style gas filtration for purifying nitrogen, and is ideal for use in LC/MS systems. The two-position baseplate (installed in the gas line) allows cartridges to be exchanged without introducing oxygen into the system. Spring-loaded check valves seal when a cartridge is removed and open only when a new cartridge has been locked in place. There is no need for loosening and tightening fittings every time you change cartridges, and your system cannot become contaminated during the changing process.

To meet the high flow needs of the LC/MS system, the charcoal-filled cartridges are positioned and connected in parallel. The incoming gas stream is split equally between the cartridges, and the two streams are rejoined after purification but before the gas exits the baseplate. This approach allows longer contact between the nitrogen and the adsorbent, ensuring higher gas purity and eliminating a potential source of contaminants to your analyses.

A handy date wheel, included with the system, indicates the cartridge installation date and the recommended replacement date.

**Table I** Super-Clean™ Filters provide high-purity outlet gas

Type of filter	Max. Flow	Outlet Gas Quality %	Maximum Pressure	Estimated Lifetime
Hydrocarbon (charcoal)	20L/min.	99.9999%	11 bar/159psi	3 to 6 months

## Super-Clean™ Gas Trapping System for LC/MS

Description	qty.	cat.#
Super-Clean™ Gas-Trapping System (2-position baseplate, 2 charcoal filters)	ea.	22062
2-Position Baseplate (1/4" Fittings)	ea.	22060
Replacement Hydrocarbon (Charcoal) Filters	2-pk.	22061



20L of purified nitrogen per minute!

GC ACCESSORIES | PURUS GAS SYSTEMS



Chris Johnston, Instrument Innovations Technician

## Super-Clean™ Click-On Traps

### Click-On Inline Super-Clean™ Traps

- High-purity output ensures 99.9999% pure gas.
- Click-On fittings for easy, leak-tight cartridge changes; brass or stainless steel, 1/4" or 1/8".
- Helium-Specific Triple Trap is ideal for GC/MS.

Using the same features and benefits as the Super-Clean™ base-plates and filters (page 220), SGT designed an inline trap. Click-On adaptor connectors allow cartridges to be exchanged without introducing oxygen. Spring-loaded check valves seal when a filter is removed and open only when a new filter has been locked in place. There is no need for loosening and tightening fittings every time a trap is changed, and your system will not become contaminated during the process.

The Triple Trap is ideal for purifying carrier gas—it contains oxygen, moisture, and hydrocarbon scrubbers in one cartridge.

The Fuel Gas Trap is ideal for purifying flame ionization detector (FID) fuel gases, removing both moisture and hydrocarbons.

The Helium-Specific Triple Trap is ideal for purifying helium in GC/MS systems. This trap is packed and purged under helium and contains oxygen, moisture, and hydrocarbon scrubbers in one cartridge.

Trap replacement depends on the quality of the incoming gas. Use the double connector and install an indicating cartridge after a trap to indicate when the trap should be replaced.



Filter Type	Gas Quality at Outlet	Maximum Pressure	Maximum Flow (L/min.)	Use For	H <sub>2</sub> O (g)	Capacity O <sub>2</sub> (mL)	Hydrocarbons (g)	Estimated Lifetime (years)
Moisture cat.#22467	>99.9999	11 bar 160psi	25	Inert carrier gas, helium, air, H <sub>2</sub>	21	NA	NA	>3
Oxygen cat.#22468	>99.9999	11 bar 160psi	25	Inert carrier gas	NA	3000	NA	>3
Hydrocarbon cat.#22466	>99.9999	11 bar 160psi	25	Inert carrier gas, helium, air, H <sub>2</sub>	NA	NA	36 <sup>3</sup>	>3
Fuel Gas <sup>1</sup> cat.#22465	>99.9999	11 bar 160psi	25	Inert carrier gas, helium, air, H <sub>2</sub>	10	NA	18 <sup>3</sup>	>2
Triple <sup>2</sup> cat.#22464	>99.9999	11 bar 160psi	25	Inert carrier gas	6	1000	12 <sup>3</sup>	>2

<sup>1</sup>Removes hydrocarbons, moisture.

<sup>2</sup>Removes hydrocarbons, moisture, oxygen.

<sup>3</sup>As *n*-butane.

**NOTE:** Super-Clean™ Gas Filters are recommended for purifying non corrosive gases with low concentrations of contaminants. The maximum concentration of oxygen in the incoming gas stream for oxygen purifiers is 0.5%.

### Click-On Inline Super-Clean™ Traps and Connector Kits

Description	qty.	cat.#
Carrier Gas Purification Kit, 1/8" Stainless Steel Includes (2) 1/8" SS connectors and (1) oxygen/moisture/hydrocarbon triple trap	kit	22456
Carrier Gas Purification Kit, 1/8" Brass Includes (2) 1/8" brass connectors and (1) oxygen/moisture/hydrocarbon triple trap	kit	22457
Carrier Gas Purification Kit, 1/4" Stainless Steel Includes (2) 1/4" SS connectors and (1) oxygen/moisture/hydrocarbon triple trap	kit	22458
Carrier Gas Purification Kit, 1/4" Brass Includes (2) 1/4" brass connectors and (1) oxygen/moisture/hydrocarbon triple trap	kit	22459
Fuel Gas Purification Kit, 1/8" Stainless Steel Includes (4) 1/8" SS connectors and (2) hydrocarbon/moisture traps	kit	22460
Fuel Gas Purification Kit, 1/8" Brass Includes (4) 1/8" brass connectors and (2) hydrocarbon/moisture traps	kit	22461
Fuel Gas Purification Kit, 1/4" Stainless Steel Includes (4) 1/4" SS connectors and (2) hydrocarbon/moisture traps	kit	22462
Fuel Gas Purification Kit, 1/4" Brass Includes (4) 1/4" brass connectors and (2) hydrocarbon/moisture traps	kit	22463



## Click-On Inline Super-Clean™ Replacement Traps

Description	qty.	cat.#
Click-On Super-Clean™ Replacement Triple Trap (removes oxygen, moisture and hydrocarbons)	ea.	22464
Click-On Super-Clean™ Replacement Fuel Gas Trap (removes moisture and hydrocarbons)	ea.	22465



## Click-On Inline Super-Clean™ Ultra-High Capacity Traps

Description	qty.	cat.#
Ultra-High Capacity Hydrocarbon Trap	ea.	22466
Ultra-High Capacity Moisture Trap	ea.	22467
Ultra-High Capacity Oxygen Trap	ea.	22468

## Helium-Specific Click-On Inline Super-Clean™ Trap and Connector Kits

Description	qty.	cat.#
<b>Kits</b> Helium-Specific Carrier Gas Cleaning Kit, 1/8" Stainless Steel Includes (2) 1/8" SS connectors and (1) oxygen/moisture/hydrocarbon helium-specific triple trap	kit	22469
Helium-Specific Carrier Gas Cleaning Kit, 1/8" Brass Includes (2) 1/8" brass connectors and (1) oxygen/moisture/hydrocarbon helium-specific triple trap	kit	22470
Helium-Specific Carrier Gas Cleaning Kit, 1/4" Stainless Steel Includes (2) 1/4" SS connectors and (1) oxygen/moisture/hydrocarbon helium-specific triple trap	kit	22471
Helium-Specific Carrier Gas Cleaning Kit, 1/4" Brass Includes (2) 1/4" brass connectors and (1) oxygen/moisture/hydrocarbon helium-specific triple trap	kit	22472
<b>Replacement Trap</b> Helium-Specific Replacement Triple Trap (removes oxygen, moisture and hydrocarbons)	ea.	22473

### did you know?

Helium-Specific Click-On Inline Super-Clean™ Trap and Kits are designed specifically for purification of helium in GC/MS systems!



## Click-On Inline Super-Clean™ Indicator

- Oxygen: Green to Grey
- Moisture: Beige to Clear

Description	qty.	cat.#
Click-On Inline Super-Clean™ Indicator (oxygen, moisture)	ea.	22474

### tech tip

Install an indicator after the Click-On inline filter so there is no confusion about when to replace the traps.



## Click-On Inline Super-Clean™ Connectors

- Click-On connectors allow you to change traps quickly, without introducing oxygen into your system.

Description	qty.	cat.#
1/8" Brass Click-On Inline Super-Clean™ Connectors	2-pk.	22475
1/8" Stainless Steel Click-On Inline Super-Clean™ Connectors	2-pk.	22476
1/4" Brass Click-On Inline Super-Clean™ Connectors	2-pk.	22477
1/4" Stainless Steel Click-On Inline Super-Clean™ Connectors	2-pk.	22478



## Click-On Inline Super-Clean™ Double Connector

- Connects any Click-On trap to a Click-On indicator.

Description	qty.	cat.#
Click-On Inline Super-Clean™ Double Connector, stainless steel	ea.	22479



## Wall-Mounting Clamps for Click-On Inline Super-Clean™ Traps

Description	qty.	cat.#
Wall-Mounting Clamps for Click-On Inline Super-Clean™ Traps	4-pk.	22480



## Replacement O-Rings for Click-On Inline Super-Clean™ Connectors

Description	qty.	cat.#
Replacement O-Rings for Click-On Inline Super-Clean™ Connectors	10-pk.	22481



## Gas-Specific Purifier Modules, Leak Detector

### VICI® Mat/Sen™ Gas-Specific Purifier Modules

- Replace separate oxygen, moisture, and hydrocarbon traps with one multiple-bed purifier, specific for purifying helium, hydrogen, nitrogen, or air.
- Reduce gas impurities from ppm to low ppb levels.
- Decrease baseline noise and increase GC/MS sensitivity.
- Pre-purged with the specified gas, to shorten downtime.



Dimensions: 21" x 1 1/2"  
(53.3 x 3.8 cm)

These gas-specific purifier modules offer dramatic reductions in contaminant levels—from ppm to levels that are below the limit of standard analytical detection. Performance is optimized by incorporating a multiple-bed format that progressively lowers concentrations of contaminants at each successive bed. VICI® Mat/Sen™ purifiers are guaranteed to produce gases that are purer than 99.9999%, when supplied with gas of 99.995% purity, and are pre-purged with the specified gas to speed conditioning. Purifier capacity is approximately four tanks of gas at 99.995% (50ppm) purity, and correspondingly longer for purer gases.

#### Specifications:

Length	21" (52.3cm)	Pressure Drop from 120psi (827 kPa) inlet at at 0-500mL/min.:	<0.20psi (2.1kPa)
Diameter	1.5" (3.8cm)	End Fittings	compression, 1/8" or 1/4", stainless steel
Maximum Inlet Pressure	1000psi (6895 kPa)	Shipping Weight	3.04 lb. (1300 g)
Maximum Recommended Flow	500mL/min.		

Please Note: We recommend using an indicating oxygen trap (e.g., cat.# 22029, pg. 220) downstream from a VICI® Mat/Sen™ purifier to continually ensure gas purity and indicate absolute change-out time.

### it's a fact

The nitrogen module is excellent for LC/MS systems.

Gas-Specific Purifier Module	Max Pressure	Compression Tube Fittings			
		1/4-inch		1/8-inch	
		qty.	cat.#	qty.	cat.#
Helium Purifier Module	1000psig, 6895 KPa	ea.	22600	ea.	22601
Hydrogen Purifier Module	1000psig, 6895 KPa	ea.	22602	ea.	22603
Nitrogen Purifier Module*	1000psig, 6895 KPa	ea.	22604	ea.	22605
Air Purifier Module	1000psig, 6895 KPa	ea.	22606	ea.	22607

\*Warning: Do not use with nitrogen containing more than 500ppm of oxygen. If the oxygen level in the stream exceeds 500ppm, use an air purifier.

### Restek Electronic Leak Detector

- Reliable thermal conductivity leak detector.
- Easy to hold and operate.
- Responds to leaks in less than 2 seconds.
- Audible alarm plus LED readout.
- Auto zeros with the touch of a button.
- Built-in rechargeable battery.



CE



Easy-to-clean probe assembly

Description	qty.	cat.#
Leak Detector with 110Volt Battery Charger	ea.	22451
Leak Detector with 220Volt European Battery Charger	ea.	22451-EUR
Leak Detector with 220Volt UK Battery Charger	ea.	22451-UK

Caution: The Restek Electronic Leak Detector is NOT designed for determining leaks of combustible gases. A combustible gas detector should be used for determining combustible gas leaks under any condition. The Restek Electronic Leak Detector may be used for determining trace amounts of hydrogen in a GC environment only.

### Leak Detector Accessory Kit

The kit includes an adaptor fitting to detect leaks in hard-to-reach locations, and a mounting bracket that can be affixed to the wall or GC.



Verify hard-to-reach leaks with the adaptor fitting.



Leak Detector Accessory Kit

Description	qty.	cat.#
Leak Detector Accessory Kit (adaptor fitting for probe, mounting bracket)	kit	22453

Leak Detector is easily accessed when stored in the mounting bracket.

**RESTEK**

**HROMalytic** **RESTEK** '07  
Australian Distributors **ECH**nology  
Tel: 03 9762 2034 Fax: 03 9761 1169 [www.chromtech.net.au](http://www.chromtech.net.au) [info@chromtech.net.au](mailto:info@chromtech.net.au)

## High-Capacity Indicating Oxygen Trap

- Indicator changes color from dark blue to black as oxygen & water are trapped.
- Lasts longer than three smaller traps.
- Use with all carrier gases.
- Ambient operating temperature, 100psig operating pressure.
- Built-in frit traps microparticles.
- Outlet gas purity:  $O_2 < 0.1\text{ppm}$  when inlet does not exceed 15ppm.  
 $H_2O < 0.5\text{ppm}$  when inlet does not exceed 10ppm.
- Includes cartridge housing and one cartridge.
- Maximum operating pressure: 150psi.
- Maximum flow: 16.5 liters/minute.



Dimensions: 9<sup>1</sup>/<sub>4</sub>" x 2"  
(23.5 x 5.1 cm)

Description	qty.	cat.#
High-Capacity Indicating Oxygen Trap, 1/8" Compression Tube Fittings	ea.	20624
High-Capacity Indicating Oxygen Trap, 1/4" Compression Tube Fittings	ea.	20623
Replacement Cartridge (fits 1/4" or 1/8" housing)	ea.	20625
Replacement O-Rings (5 small O-rings and 5 large O-rings)	kit	22081



## Indicating Oxygen Trap

- Indicator changes color from light green to grey as oxygen is trapped.
- Heavy-walled glass body prevents oxygen & water infusion.
- Pre-purged for fast stabilization.
- 100psig maximum operating pressure.
- Reduces oxygen to 0.1ppm.

Description	qty.	cat.#
Indicating Oxygen Trap, 1/8" Fittings	ea.	22010
Indicating Oxygen Trap, 1/4" Fittings	ea.	22011



Dimensions: 10" x 1<sup>1</sup>/<sub>4</sub>"  
(25.4 x 3.2 cm)

## High-Capacity Oxygen Trap

- Removes up to 600mg of oxygen or 2g of water.
- Long life—typically purifies more than five 200ft.<sup>3</sup> cylinders.
- Reduces oxygen to 15ppb.
- Maximum operating pressure: 250psi.
- Flow: 3 liters/minute @ 32psi.

Description	qty.	cat.#
High-Capacity Oxygen Trap, 1/8" Fittings	ea.	20601
High-Capacity Oxygen Trap, 1/4" Fittings	ea.	20600



Dimensions: 11" x 1<sup>1</sup>/<sub>2</sub>"  
(27.9 x 3.8 cm)

## Rechargeable Molecular Sieve S-Trap

- Traps water vapor; increases column and oxygen trap lifetime.
- Reduces baseline noise from sensitive detectors such as ECDs and mass spectrometers.
- Activated and ready to use.
- Reduces water to less than 1ppm.
- Fits in GC oven for easy thermal recharging.
- Maximum flow: 1 liter/minute.

Description	qty.	cat.#
Rechargeable Molecular Sieve S-Trap, 1/8" Fittings	ea.	20686
Rechargeable Molecular Sieve S-Trap, 1/4" Fittings	ea.	20685



Dimensions: 6<sup>3</sup>/<sub>4</sub>" x 5<sup>5</sup>/<sub>8</sub>"  
(17.1 x 14.3 cm)

## tech tip

### Carrier Gas Purity

Carrier gas should contain less than 1ppm of oxygen, moisture, or other trace contaminants, to prevent column degradation, increase column lifetime, and decrease stationary phase bleed.

The expense of using high-purity gases in combination with carrier gas line purifiers will be offset by longer column lifetime and less GC maintenance.

## Moisture Traps, Hydrocarbon Traps



Dimensions: 11" x 1 1/2"  
(27.9 x 3.8 cm)

### High-Capacity Moisture Trap

- Purged with ultra-high-purity helium; ready to use.
- Reduces water to less than 15ppb.
- Maximum operating pressure: 250psi.
- Maximum flow: 1.25 liters/minute.

Description	qty.	cat.#
High-Capacity Moisture Trap, 1/8" Fittings	ea.	21997
High-Capacity Moisture Trap, 1/4" Fittings	ea.	20638



Dimensions: 13" x 2"  
(33 x 5.1 cm)

### Indicating Moisture Trap

- Reduces water to less than 10ppb; indicator changes color from blue to pink at 5% relative humidity.
- Pre-purged for fast stabilization.
- Reduces noise from high-sensitivity detectors.
- Heavy-walled glass body prevents oxygen & water infusion.
- 40µm frit prevents microparticulate damage to needle valves and flow controllers.
- Maximum operating pressure: 100psig.

Description	qty.	cat.#
Indicating Moisture Trap, 1/8" Fittings	ea.	22014
Indicating Moisture Trap, 1/4" Fittings	ea.	22015



Dimensions: 11" x 1 1/2"  
(27.9 x 3.8 cm)

### Capillary-Grade Hydrocarbon Trap

- Packed with an extremely high surface area, baked coconut shell-based activated carbon.
- Purged with ultra-high-purity helium.
- Reduces organics to 0.1ppm (assuming 100ppm input).
- Maximum operating pressure: 250psi.

Description	qty.	cat.#
Capillary-Grade Hydrocarbon Trap, 1/8" Fittings	ea.	21991
Capillary-Grade Hydrocarbon Trap, 1/4" Fittings	ea.	21992



Dimensions: 9 1/4" x 2 1/4"  
(23.5 x 5.7 cm)

### Refillable Hydrocarbon Trap

- Removes trace impurities from carrier gas: reduces organics to 0.1ppm (assuming 100ppm input).
- 20µm frit prevents gas contamination by purifier particles.
- Stops carrier gas contaminants from interfering with purge & trap systems.
- Refillable and rechargeable.
- Maximum operating pressure: 250psi.
- Maximum flow: 1.25 liters/minute.

Description	qty.	cat.#
Refillable Hydrocarbon Trap, 1/8" Fittings	ea.	22012
Refillable Hydrocarbon Trap, 1/4" Fittings	ea.	22013
Carbon Refill (two recharges)	pint	20626



Dimensions: 6 3/4" x 5 5/8"  
(17.1 x 14.3 cm)

### Hydrocarbon S-Trap

- Removes hydrocarbons and other contaminants from the gas stream.
- Reduces organics to 0.1ppm (assuming 100ppm input).
- Each trap individually activated to ensure maximum efficiency.
- Fits in GC oven for easy thermal recharging.
- Maximum operating pressure: 60psi.

Description	qty.	cat.#
Hydrocarbon S-Trap, 1/8" Fittings	ea.	22016



Dimensions: 6" x 1 3/4"  
(15.2 x 4.4 cm)

### Indicating Hydrocarbon Trap for Air Compressors

- Pass compressed air from an oil-filled air compressor through this trap, to remove oil vapors and mist.
- Indicator changes from pale pink to deep pink.

Description	qty.	cat.#
Indicating Hydrocarbon Trap for Air Compressors, 1/8" Fittings	ea.	20637
Indicating Hydrocarbon Trap for Air Compressors, 1/4" Fittings	ea.	20636

## High-Capacity Split Vent Trap

- Reduces the release of hazardous materials from the capillary split vent into the lab.
- Lasts approximately one month or 1,500 injections.
- Includes connecting lines and mounting kit.

Description	qty.	cat.#
High-Capacity Split Vent Trap	ea.	20698
High-Capacity Split Vent Trap	5-pk.	20699



Dimensions: 6" x 1"  
(15.2 x 2.5 cm)

## ECD Vent Trap

- Reduces the release of hazardous materials from the ECD vent into the lab.
- Includes connecting lines and mounting kit.

Description	qty.	cat.#
ECD Vent Trap	ea.	22017



Dimensions: 6" x 1"  
(15.2 x 2.5 cm)

## Replacement Chemical Traps for Agilent GCs


- Easy to install.
- Attach to same fittings as original equipment.
- Built-in frits retain fine particles; adsorbents remove both moisture and hydrocarbons.

Description	Similar to Agilent part #	qty.	cat.#
A) Replacement Chemical Trap for Agilent 6890/6850 GCs	G1544-80550	ea.	22820
B) Replacement Chemical Trap for Agilent 5890 GCs	05890-61260	ea.	21610
C) Replacement Chemical Trap for Agilent 5880 GCs	19362-60500	ea.	21158
D) Split Vent Line for Agilent GCs (32-inch) (Includes all installation hardware.)	19251-80525 G1544-20620	2-pk.	22800
E) O-rings for Agilent Trap Fittings	5180-4181	25-pk.	22064
F) Optional Split Vent Trap Assembly for Agilent 6890/6850 GCs	G1544-60610	kit	23031
G) Replacement Traps (2) and O-rings (4)	G1544-80530	kit	23032




## tech tip

Eliminate ghost peaks—change your chemical trap often!



## tech tip

Gas generators are an economical source of pure gases, and eliminate the inconvenience and danger of high pressure cylinders.



## a plus 1 story

"We enjoy doing business with Restek. Their technical knowledge and willingness to back their products help us to maximize the performance of our chromatography instruments."

**Jean-François Vergelin,**  
Département de Seine et Marne, Direction de l'Eau et de l'Environnement, Laboratoire Départemental d'Analyse des eaux (Melun, France)

## Hydrogen Generators



### Parker Gas Generators

More labs than ever are switching from high-pressure cylinders to generated gas. In evaluating the shift to on-site gas generation, you should consider convenience, reliability, safety, and cost as your primary criteria. Restek has developed a marketing alliance with Parker Hannifin Corporation, Filtration and Separation Division—formerly Whatman—to provide the best value in laboratory gas generators worldwide. Produced and supported by an ISO 9001 registered organization, Parker Balston's hydrogen generators are built to meet the toughest laboratory standards—CSA, UL, CE, and IEC 1010.



### Parker Hydrogen Generators

- Selectable delivery pressure: 0–100psig.
- High hydrogen purity—99.9995%—for better chromatography.
- No high-pressure cylinders—greater convenience and improved lab safety.

Fuel-grade purity hydrogen generators are safer alternatives to high-pressure gas cylinders. An exclusive solid polymer electrolyte produces hydrogen on demand. Deionized water and an electrical supply are all you need to generate hydrogen for weeks of continuous operation. With an output capacity of up to 500cc/min., one generator can supply 99.9995% pure carrier gas to several chromatographs and fuel gas to 12 FIDs. Based on cylinder gas savings alone, a hydrogen generator will pay for itself in less than one year.

Parker hydrogen generators are certified for laboratory use by Canadian Standards Association (CSA), Underwriters Laboratories (UL), and International Electrotechnical Commission (IEC) 1010. A built-in sensing circuit shuts the generator down if a hydrogen leak is detected.

#### Specifications

Purity:	99.9995% pure hydrogen
	oxygen < 1ppm
	moisture < 1ppm
	hydrocarbons/halocarbons: 1ppb
Delivery Pressure:	2 to 30psig $\pm$ 3%
	30 to 100psig $\pm$ 2%
Pressure Control:	5 to 20psig $\pm$ 0.5%
	20 to 100psig $\pm$ 0.2%
Connections:	1/8" tube
Electrical Requirements:	117 VAC/234 VAC
Physical Dimensions:	14.75" h x 13" w x 14" d
	(37.8 x 33.3 x 35.9cm)
Shipping Weight:	40 lbs. (18kg)

### free literature

For more information on the benefits of using hydrogen generators, download a free copy of our *ChromGas® Hydrogen Generators* from Parker Fast Facts from [www.restek.com](http://www.restek.com)!

lit. cat.# 580053

### ordering note

For **international orders**, please add the appropriate power cord suffix from the table below.

Description	Capacity	qty.	cat.#
Hydrogen Generator A9090	90cc/min.	ea.	22033
Hydrogen Generator A9090 with European Power Cord	90cc/min.	ea.	22033-551
Hydrogen Generator A9150	160cc/min.	ea.	22034
Hydrogen Generator A9150 with United Kingdom Power Cord	160cc/min.	ea.	22034-550
Hydrogen Generator B9200	250cc/min.	ea.	22035
Hydrogen Generator B9400	500cc/min.	ea.	22036
<b>Replacement Components for Hydrogen Generators (for all models listed above)</b>			
Replacement Deionizer Bag		2-pk.	21670
Replacement Desiccant Cartridge		ea.	21671

### International Power Cord Sets

Just add the proper suffix to the catalog number for the gas generator you are ordering.

Location	qty.	cat.# suffix
United Kingdom (230VAC, 50/50Hz)	ea.	-550
European (230VAC, 50/60Hz)	ea.	-551
IEC Connector Only (230VAC, 50/60Hz)	ea.	-552
Japanese (200VAC, 50/60Hz)	ea.	-556
Japanese for Zero Air (100VAC, 50/60Hz)	ea.	-553
Japanese for Hydrogen (100VAC, 50/60Hz)	ea.	-554
Japanese for Nitrogen (100VAC, 50/60Hz)	ea.	-555



**Model FID-1000 and FID-2500 Gas Stations**

- Single unit produces UHP zero air from house compressed air and 99.9995% pure hydrogen from deionized water.
- Ideal for supplying up to 5 - 6 FIDs.
- Eliminates inconvenient and dangerous cylinders.
- Silent operation, minimal operator attention required.

Parker Balston's Gas Stations provide both UHP grade hydrogen fuel gas and zero grade air for flame ionization detectors on gas chromatographs. The system is specifically designed to supply fuel gas to FIDs and to support flame thermionic and flame photometric detectors. The units produce zero air by purifying compressed air to a total hydrocarbon concentration of 0.1ppm or less (measured as methane).

The hydrogen generators produce hydrogen gas from deionized water, using the principle of electrolytic dissociation of water and hydrogen proton conduction through a proton exchange membrane cell.

**Specifications**

Hydrogen Purity:	99.9995%
Zero Air Purity:	FID-1000: < 0.1ppm total hydrocarbons as methane FID-2500: < 0.05ppm total hydrocarbons as methane
Max. Hydrogen Flow Rate:	FID-1000: 90cc/min. FID-2500: 250cc/min.
Max. Zero Air Flow Rate:	FID-1000: 1000cc/min. FID-2500: 2500cc/min.
Power:	120VAC/amp, 60Hz, 400 watts
Hydrogen Outlet Pressure:	60 psig
Zero Air Outlet Pressure:	40-125 psig*
Inlet Connection:	1/4" NPT (female)
Outlet:	1/8" compression
Dimensions:	16.5"h x 10.5"w x 17"d (42cm x 27cm x 43cm)
Weight:	53 lbs. (24kg)



Produce zero air and pure hydrogen from one unit!

GC ACCESSORIES | PURUS GAS SYSTEMS


**free literature**

For more information, download a free copy of our *FID Gas Stations* Fast Facts from [www.restek.com](http://www.restek.com)!

lit. cat.# 580051

**ordering note**

For **international orders**, please add the appropriate power cord suffix from the table below.



Description	qty.	cat. #
Model FID-1000 Gas Station (ideal for 1 - 2 FIDs)	ea.	20177
Model FID-2500 Gas Station (ideal for 5 - 6 FIDs)	ea.	24913
<b>Replacement Components for FID Gas Stations</b>		
Resin Bed Cartridge for Hydrogen Generators in FID-1000 and FID-2500 Gas Stations	ea.	24914
Replacement Desiccant Cartridge	ea.	21671
FID Gas Station Maintenance Kit (Includes 1 desiccant cartridge, 1 resin bed cartridge, 1 filter cartridge)	ea.	24915

\*Zero air inlet requires minimum of 40psig compressed air pressure.

**International Power Cord Sets**

Just add the proper suffix to the catalog number for the gas generator you are ordering.

Location	qty.	cat.# suffix
United Kingdom (230VAC, 50/50Hz)	ea.	-550
European (230VAC, 50/60Hz)	ea.	-551
IEC Connector Only (230VAC, 50/60Hz)	ea.	-552
Japanese (200VAC, 50/60Hz)	ea.	-556
Japanese for Zero Air (100VAC, 50/60Hz)	ea.	-553
Japanese for Hydrogen (100VAC, 50/60Hz)	ea.	-554
Japanese for Nitrogen (100VAC, 50/60Hz)	ea.	-555

## Zero Air Generators



### Parker Balston® Zero Air Generators

- Turn in-house compressed air into ultra-pure air (<0.1ppm total hydrocarbons).
- Remove hydrocarbons to less than 0.1ppm by catalytic oxidation.
- Operate at 40 to 125psi.
- Typical payback is less than one year, based on cylinder costs.
- Install easily and take up little bench space.
- Maintenance kits include a one year supply of prefilters and final filter.

#### Specifications

Maximum Zero Air Flow Rate:	75-83NA	1 lpm
	HPZA-3500	3.5 lpm
	HPZA-7000	7 lpm
	HPZA-18000	18 lpm
	HPZA-30000	30 lpm
Outlet Hydrocarbon Concentration (as methane):	75-83NA	< 0.1 ppm
	HPZA-30000	< 0.1 ppm
	Other Models	< .05 ppm
Minimum/Maximum Inlet Air Pressure:	40 psig / 125 psig	
Maximum Inlet Hydrocarbon Concentration (as methane):	100 ppm	
Pressure Drop at Maximum Flow Rate:	4 psid	
Maximum Inlet Air Temperature:	78°F (25°C)	
Inlet/Outlet Ports:	1/4" NPT (female)	
Start-up Time to Specified Hydrocarbon Concentration:	45 minutes	
Electrical Requirements:	75-83NA	120 VAC/60 Hz, 0.5 amps
	Other Models	120 VAC/60 Hz, 3.5 amps
Dimensions:	75-83NA	12"h x 10"w x 3"d (30cm x 25cm x 8cm)
	Other Models	16"h x 11"w x 13"d (42cm x 27cm x 34cm)
Shipping Weight:	75-83NA	7 lbs. (3 kg)
	Other Models	41 lbs. (19 kg)

Model	Number of FIDs*
75-83NA	Up to 3
HPZA-3500	Up to 11
HPZA-7000	Up to 23
HPZA-18000	Up to 60
HPZA-30000	Up to 100

\*based on a 300 cc/min. fuel air rate

Zero Air Generator	Capacity	qty.	cat. #
Zero Air Generator Model 75-83NA	1000cc/min.	ea.	20684
Zero Air Generator Model 75-83NA with United Kingdom Power Cord	1000cc/min.	ea.	20684-550
Zero Air Generator Model HPZA-3500	3500cc/min.	ea.	20680
Zero Air Generator Model HPZA-3500 with European Power Cord	3500cc/min.	ea.	20680-551
Zero Air Generator Model HPZA-7000	7000cc/min.	ea.	20681
Zero Air Generator Model HPZA-18000	18,000cc/min.	ea.	20682
Zero Air Generator Model HPZA-30000	30,000cc/min.	ea.	20683
<b>Maintenance Kits (includes a one-year supply of prefilters and final filter)</b>		<b>qty.</b>	<b>cat. #</b>
Maintenance Kit for Model 75-83NA		kit	21646
Maintenance Kit for Models HPZA-3500, HPZA-7000, HPZA-18000, HPZA-30000		kit	21647
<b>Replacement Catalyst Towers</b>	<b>Capacity</b>	<b>qty.</b>	<b>cat. #</b>
Replacement Catalyst Tower for Model 75-83NA	1000cc/min.	ea.	22005
Replacement Catalyst Tower for Model HPZA-3500	3500cc/min.	ea.	22004
Replacement Catalyst Tower for Model HPZA-7000	7000cc/min.	ea.	22006
Replacement Catalyst Tower for Model HPZA-18000	18,000cc/min.	ea.	22007
Replacement Catalyst Tower for Model HPZA-30000	30,000cc/min.	ea.	22008

### ordering note

For **international orders**, please add the appropriate power cord suffix from the table below.

### free literature

For more information, download a free copy of our *Zero Air Generators Fast Facts* from [www.restek.com](http://www.restek.com)!

lit. cat.# 580050

### International Power Cord Sets

Just add the proper suffix to the catalog number for the gas generator you are ordering.

Location	qty.	cat.# suffix
United Kingdom (230VAC, 50/50Hz)	ea.	-550
European (230VAC, 50/60Hz)	ea.	-551
IEC Connector Only (230VAC, 50/60Hz)	ea.	-552
Japanese (200VAC, 50/60Hz)	ea.	-556
Japanese for Zero Air (100VAC, 50/60Hz)	ea.	-553
Japanese for Hydrogen (100VAC, 50/60Hz)	ea.	-554
Japanese for Nitrogen (100VAC, 50/60Hz)	ea.	-555



## Parker Nitrogen Generators

- Turn compressed air into ultra-pure nitrogen (up to 99.9995%).
- Flows from 1 to 75+ lpm.
- Require only a compressed air source and 110 volt AC power.
- Safe, reliable, low maintenance.
- Maintenance kits include replacement filters.

### Specifications

	Model HPN2-1100 or UHPN2-1100	Model HPN2-2000	Model N2-14 or N2-14A
Maximum Nitrogen Flow Rate:	See Flow Table	2 lpm	78scfh** at 95% purity
Nitrogen Purity:	99.9999%	99.99%	95.0%–99.5%
Maximum Nitrogen Outlet Pressure:	See Flow Table	90 psig	
CO Concentration:	< 1.0 ppm	NA	
CO <sub>2</sub> Concentration:	< 1 ppm	< 1 ppm	
O <sub>2</sub> Concentration:	< 1 ppm	< 100 ppm	
H <sub>2</sub> O Concentration:	≤ 2 ppm	≤ 2 ppm	
Hydrocarbon Concentration <sup>1</sup> :	< 0.1 ppm	NA	
Argon Concentration <sup>2</sup> :	0.9%	0.9%	
Atmospheric Dewpoint:			-58°F (-50°C)
Suspended Liquids:			None
Particles > 0.01µm:			None
Oxygen Analyzer:			Included with Model 75-720NA
Commercially Sterile:			Yes
Minimum/Maximum Inlet Pressure:	60 psig/125 psig	75 psig/120 psig	60/145 psig
Maximum Pressure Drop (99% N <sub>2</sub> Purity, 125 psig):			10 psig
Recommended Inlet Temperature:	≤ 78°F (25°C)	≤ 78°F (25°C)	≤ 68°F (25°C) (Max.)
Ambient Operating Temperature:	60°F–100°F (16°C–38°C)	40°F–100°F (4°C–38°C)	110°F (43°C) (Max.)
Maximum Air Consumption:	42 lpm (1.5 scfm)*	42 lpm (1.5 scfm)*	
Inlet Connection:	1/4" NPT (female)	1/4" NPT (female)	1/4" NPT
Outlet Connection:	1/4" compression	1/8" NPT compression	1/8" NPT
Electrical Requirements <sup>3</sup> :	120 VAC/60 Hz	120 VAC/60 Hz	N2-14: None N2-14A: 120 VAC/60 Hz/25 Watts
Dimensions:	35"h x 12"w x 16"d (89cm x 30cm x 41cm)	35"h x 12"w x 16"d (89cm x 30cm x 41cm)	50"h x 16"w x 16"d (127cm x 41cm x 41cm)
Shipping Weight:	115 lbs. (52 kg)	115 lbs. (52 kg)	N2-14: 75 lbs. (34 kg) N2-14A: 80 lbs. (36 kg)



<sup>1</sup>Models HPN2-1100 and HPN2-2000 do not remove hydrocarbons.  
<sup>2</sup>Purity specification for nitrogen does not include argon concentration.  
<sup>3</sup>Power consumption is:  
 Model HPN2-1100 = 25 Watts  
 Model UHPN2-1100 = 700 Watts  
 Model HPN2-2000 = 25 Watts

### Flow Table for Models HPN2-2000, HPN2-1100, and UHPN2-1100

Inlet Air Pressure (psig)	Maximum Outlet Flow (cc/min.)	Maximum Outlet Pressure (psig)
	<i>Models HPN2-1100 and UHPN2-1100</i>	
125	1100	85
110	1000	75
100	900	65
90	800	60
80	700	50
70	600	45
60	500	35
	<i>Model HPN2-2000</i>	
75-120	2000	90

### Nitrogen Generator

	qty.	cat.#
Nitrogen Generator N2-14 (general purpose) 78 scfh** max. flow at 95% purity	ea.	20677
Nitrogen Generator N2-14 with European Power Cord	ea.	20677-551
Nitrogen Generator N2-14A (general purpose w/oxygen analyzer) 78 scfh** max. flow at 95% purity	ea.	21652
Nitrogen Generator HPN2-2000 (high purity—99.99%) 2.0 lpm max. flow	ea.	21654
Nitrogen Generator HPN2-1100 (ultra-high purity—99.9995%) 1.1 lpm max. flow	ea.	21653
Nitrogen Generator UHPN2-1100 (ultra-high purity—99.9995%); [HC<0.1ppm] 1.1 lpm max. flow	ea.	20697

### Maintenance Kits

	qty.	cat.#
Maintenance Kit for Models N2-14, N2-14A, 75-72, 75-720NA	kit	21648
Maintenance Kit for Models HPN2-1100, HPN2-2000, 76-96, 76-92	kit	21649
Maintenance Kit for Models UHPN2-1100, 76-94	kit	21655

\*Standard cubic feet per minute.

\*\*Standard cubic feet per hour.

## free literature

For more information, download a free copy of our *Nitrogen Generators Fast Facts* from [www.restek.com](http://www.restek.com)!

lit. cat.# 580052

## ordering note

For **international orders**, please add the appropriate power cord suffix from the table on the previous page.

## Pressure Regulators: Introduction

### Introduction to Pressure Regulators



#### General Purpose or Analytical?

General-purpose regulators usually are best suited for applications involving gases that are less than 99.995% pure: pneumatically-actuated valves and autosamplers, blanketing, inert atmospheres, and any other application not directly integrated with analytical data production. General purpose regulators have nylon-reinforced neoprene diaphragms that provide very good pressure control but are prone to air and moisture diffusion and hydrocarbon off-gassing.

Analytical regulators are recommended for applications in which maintaining the purity of a gas or mixture is the overriding concern, i.e., for applications requiring gases that are greater than 99.995% pure. They are commonly used in analytical labs. Analytical regulators have stainless steel diaphragms for pressure control. Stainless steel is not subject to the diffusion and off-gassing associated with neoprene diaphragms, and is easily purged of atmospheric contaminants when put into service.

#### Dual- or Single-Stage?

Dual-stage regulators reduce the source pressure to outlet pressure in two steps. The first stage reduces the inlet pressure to about three times the maximum working pressure. Outlet pressure regulation is controlled by the second stage and is set through the use of an adjusting knob. This two-step regulation is highly recommended for services requiring a near constant delivery pressure as the source pressure decays, including chromatographic analyses.

Single-stage regulators perform the same function as dual-stage regulators, but in a single step-down from source pressure to outlet pressure. For this reason, the outlet pressure cannot be as accurately maintained as the source pressure decays. We highly recommend that single-stage regulators be used only in circumstances in which the operator can monitor and adjust the regulator as needed, when the regulator is supplied with a nearly constant source pressure, or when additional pressure regulation is supplied downstream.

#### Brass or Stainless Steel?

Analytical regulators made from brass bar stock provide optimum performance for most analytical applications. Brass provides excellent strength and cleanliness and the machined bar stock design has less dead volume than forged-body regulators, making purging of atmospheric contaminants faster and more assured.

Regulators with stainless steel bodies were designed for delivering corrosive gases that would be incompatible with brass. With the advent of semiconductor manufacturing and high sensitivity analytical techniques, stainless steel also has proven to be a better surface for removing “sticky” atmospheric contaminants that interfere with detectors downstream. Unless these regulators are used in an all-stainless-steel system that incorporates welded tubing and special fittings, and in which rigorous cleaning and proper gas management are practiced, the extra expense relative to brass is not justified.

## Overview of Restek's Brass and Stainless Steel Body Ultra-High-Purity Regulators

These regulators feature metal-to-metal seals throughout for long-term leak-tightness, and a metal diaphragm outlet valve ensures gas purity. Each regulator is helium leak-test-certifiable to  $1 \times 10^{-8}$  scc/sec. and is fully assembled and tested for your convenience. 100psig maximum delivery pressure supports EPC operation. Maximum inlet pressure is 3000psig. Brass bar stock construction minimizes dead volume. Stainless steel construction is more easily purged of atmospheric contaminants, and is more resistant to attack from dry corrosive gases.

### Ultra-High-Purity Brass Body Regulators

These regulators are the best choice when using ultra-high-purity carrier gas for sensitive GC applications using MS, PID, HID, or ECD detection methods. They feature reduced internal dead-volume, relative to stainless steel bodies. The stainless steel regulator diaphragm minimizes non-metal parts, to help ensure ultra-high gas purity, and the metal valve diaphragm ensures leak-free shut-off. Oxidation-resistant chrome plating maintains a like-new appearance.

### Dual-Stage Ultra-High-Purity Brass Regulators

- Oxidation-resistant, chrome-plated.
- Most stable outlet pressure control throughout the life of a high-pressure gas cylinder.
- Secondary pressure regulation not needed.
- Most widely used regulator.
- Less internal volume than stainless steel regulators.

Outlet pressure: 0 to 100psig  
 Outlet gauge: 30" – 0 to 200psig  
 Inlet gauge: 0 to 4000psig  
 Outlet assembly: diaphragm valve, 1/4" tube fitting

Fitting	qty.	cat.#
CGA 580 (N <sub>2</sub> , He, Ar)	ea.	21667
CGA 350 (H <sub>2</sub> , P <sub>2</sub> )	ea.	21668
CGA 590 (Air)	ea.	21669



### Single-Stage Ultra-High-Purity Brass Regulators

- Oxidation-resistant, chrome-plated.
- Use when there is secondary pressure regulation downstream.
- Identical gas purity protection as with our dual-stage regulators.

Outlet pressure: 0 to 100psig  
 Outlet gauge: 30" – 0 to 200psig  
 Inlet gauge: 0 to 4000psig  
 Outlet assembly: diaphragm valve, 1/4" tube fitting

Fitting	qty.	cat.#
CGA 580 (N <sub>2</sub> , He, Ar)	ea.	20646
CGA 350 (H <sub>2</sub> , P <sub>2</sub> )	ea.	20647
CGA 590 (Air)	ea.	20648



### Ultra-High-Purity Brass Line Regulator

- Oxidation-resistant, chrome-plated.
- Use where you need to reduce the line pressure by 20psi or more.
- Same purity protection as high-pressure cylinder regulators.

Inlet connections: 1/4" FPT  
 Outlet assembly: 1/4" FPT port

Fitting	Outlet Gauge	Outlet Pressure	qty.	cat.#
1/4" female NPT ports*	30" - 0 to 100psig	0-50psig	ea.	21666
1/4" female NPT ports*	30" - 0 to 200psig	0-100psig	ea.	22452

\*Order appropriate male connector, pipe-to-tube fittings.



### Male Connector, Pipe-to-Tube Fittings

Fitting Type	Size (inches)	Parker #	Similar to Swagelok®	Brass		Stainless Steel	
				qty.	cat.#	qty.	cat.#
Male Connector	1/4" to 1/4" NPT	4 MSC 4N	400-1-4	10-pk.	21842	2-pk.	21942
Male Connector	1/8" to 1/4" NPT	2 MSC 4N	200-1-4	10-pk.	21844	2-pk.	21944
Tube End Reducer	1/4" tube to 1/8"	4 TUR 2	200-R-4	5-pk.	21834	2-pk.	21934



male connector



tube end reducer

## Stainless Steel Pressure Regulators, Switchover Systems

### Ultra-High-Purity Stainless Steel Body Regulators

These regulators are the standard for ultra-high-purity and corrosion-resistant pressure regulation. They are more easily purged of atmospheric components, compared to brass regulators, making them ideal for the most demanding applications. Regulation performance is equal to our brass body regulators. For use in all-stainless steel systems where welded tubing and special fittings are used, and rigorous cleaning and proper gas management are practiced.

### Dual-Stage Ultra-High-Purity Stainless Steel Regulators

- Most stable outlet pressure control throughout the life of a high-pressure gas cylinder.
- Secondary pressure regulation not needed.

Outlet pressure: 0 to 100psig  
 Outlet gauge: 30" – 0 to 200psig  
 Inlet gauge: 0 to 4000psig  
 Outlet assembly: diaphragm valve, 1/4" tube fitting

Fitting	qty.	cat.#
CGA 580 (N <sub>2</sub> , He, Ar)	ea.	20662
CGA 350 (H <sub>2</sub> , P <sub>2</sub> )	ea.	20663
CGA 590 (Air)	ea.	20664

### Single-Stage Ultra-High-Purity Stainless Steel Regulators

- Use when there is secondary pressure regulation downstream.
- Identical gas purity protection as with our dual-stage regulators.

Outlet pressure: 0 to 100psig  
 Outlet gauge: 30" – 0 to 200psig  
 Inlet gauge: 0 to 4000psig  
 Outlet assembly: diaphragm valve, 1/4" tube fitting

Fitting	qty.	cat.#
CGA 580 (N <sub>2</sub> , He, Ar)	ea.	20665
CGA 350 (H <sub>2</sub> , P <sub>2</sub> )	ea.	20666
CGA 590 (Air)	ea.	20667

### Critical Purity Automatic Switchover System for Non-Corrosive Service

High-purity automatic switchover systems provide a continuous supply of high purity gas to the laboratory, process, or instrument, to allow you to replace a depleted gas source without interruption in the supply of gas. Continuous supply is achieved by setting the two regulators at slightly different pressures, to discharge one side of the system at a time. These models include flexible, all-stainless-steel pigtailed with armor casing. The CGA connection on each pigtail has a check valve in the gland to prevent contamination and minimize purging requirements.

Brass Automatic Switchover System with Line Regulator	qty.	cat.#
CGA 580 (N <sub>2</sub> , He, Ar)	ea.	20668580
CGA 350 (H <sub>2</sub> , P <sub>2</sub> )	ea.	20668350
CGA 590 (Air)	ea.	20668590
Stainless Steel Automatic Switchover System with Line Regulator	qty.	cat.#
CGA 580 (N <sub>2</sub> , He, Ar)	ea.	21593580
CGA 350 (H <sub>2</sub> , P <sub>2</sub> )	ea.	21593350
CGA 590 (Air)	ea.	21593590

### Protocol Station

The protocol station is designed for convenient wall mounting of high-purity regulators. Wall mounting provides ease of use, prevents regulator damage, and improves safety. Either chrome-plated brass or 316 stainless steel option is complete with a 3-foot, flexible, all-stainless-steel pigtail with armor casing. The CGA connection on the pigtail has an integral check valve in the gland to prevent contamination during cylinder changeout.

Chrome-Plated Brass Protocol Station*	qty.	cat.#
CGA 580 (N <sub>2</sub> , He, Ar)	ea.	21347
CGA 350 (H <sub>2</sub> , P <sub>2</sub> )	ea.	21348
CGA 590 (Air)	ea.	21349
Stainless Steel Protocol Station*	qty.	cat.#
CGA 580 (N <sub>2</sub> , He, Ar)	ea.	21327
CGA 350 (H <sub>2</sub> , P <sub>2</sub> )	ea.	21328
CGA 590 (Air)	ea.	21329

\*Pressure regulator not included. Order separately.



**Switching pressure:**  
200psig/170psig  
**Inlet connections:**  
flexible SS pigtailed (36")  
**Line regulator:** 0 to 100psig



**CGA Fittings**

CGA-specified nuts and nipples with internal frit, 1/4-inch NPT nickel-plated brass.

Description	qty.	cat.#
CGA 580 Fitting, (N <sub>2</sub> , He, Ar)	ea.	21336
CGA 350 Fitting, (H <sub>2</sub> , P <sub>2</sub> )	ea.	21337
CGA 590 Fitting, (Air)	ea.	21338



**Flexible Stainless Steel Hoses**

Description	qty.	cat.#
Flexible Stainless Steel Hose, 36", 1/4" Female NPT Fittings	ea.	21339
Flexible Stainless Steel Hose, 18", 1/4" Female NPT Fittings	ea.	21340
Flexible SS Hose, 36", with Stainless Steel CGA 580	ea.	21344
Flexible SS Hose, 36", with Stainless Steel CGA 350	ea.	21345
Flexible SS Hose, 36", with Stainless Steel CGA 590	ea.	21346



**Flammable Gas Flash Arrestor—Factory Mutual Approved\***

- Gas flow shuts off in the event of flashback.
- Flame extinguished—flame front prevented from reaching the gas supply.
- No gas flow restriction under normal operating conditions.

Description	qty.	cat.#
Flammable Gas Flash Arrestor, Brass Body	ea.	21334

\*Approved for brass body servicing hydrogen, acetylene, propane, or natural gas only.



**Backpressure Regulator**

Capillary GC inlet systems have backpressure regulators to maintain a constant upstream pressure and rapidly respond to catastrophic leaks. The 0–60psig operating range is sufficient to operate a 105m, 0.25mm ID column at its optimum flow rate.

Description	qty.	cat.#
Backpressure Regulator	ea.	20635



**MINICYL Regulator**

This compact general purpose regulator has many laboratory applications including air-drying glassware, sparging or evaporating solutions, and controlling pneumatic valves. It is constructed of lightweight aluminum with an elastomer diaphragm. Includes a 0–60psig gauge and either 1/8- or 1/4-inch tube fittings.

Description	qty.	cat.#
MINICYL Regulator, 1/8" Fittings	ea.	20610
MINICYL Regulator, 1/4" Fittings	ea.	20611



**Cylinder Valve Wrench**

This specially-designed wrench enables easy opening of cylinder valves that are fitted with a hand wheel. It is also suitable for removing difficult cylinder caps.

Description	qty.	cat.#
Cylinder Valve Wrench	ea.	21321



**Universal Cylinder Wrench**

Use this versatile wrench for tightening gauges and regulator CGA fittings to cylinder outlets and pipe thread connections.

Description	qty.	cat.#
Universal Cylinder Wrench	ea.	21322



**Wall-Mount Cylinder Holder**

Don't risk injury or damage—this holder is designed to prevent free-standing cylinders from tipping over. The cast aluminum holder body can be secured to a wall or the side of a work bench. Will secure a cylinder from 4–15 inches in diameter.

Description	qty.	cat.#
Wall-Mount Cylinder Holder	ea.	21333



## Instrument-Grade Tubing



### ordering note

Required length in meters  
x 3.2808 = length in feet.

### Instrument-Grade Welded and Drawn 304 Grade Stainless Steel Tubing

Clean tubing is critical to the delivery of pure gas to your instrument. Restek's Stainless Steel Tubing is specially cleaned for inertness by using the procedure for processing our Silcosteel® and Siltek®-treated products, because scrupulously clean parts are a prerequisite for a quality coating.

Tubing Dimensions		Length (per-foot pricing on 26 feet or more)		
ID (in.)	OD (in.)	25 Feet	26–100 Feet	>100 Feet*
0.01" ID	3/16" OD	cat.# 21500	cat.# 21501	cat.# 21502
0.02" ID	1/16" OD	21503	21504	21505
0.03" ID	3/16" OD	21506	21507	21508
0.04" ID	1/16" OD	21509	21510	21511
0.085" ID	3/8" OD	21512	21513	21514
0.21" ID	1/4" OD	21515	21516	21517

\*The availability of long lengths (continuous lengths up to 500 feet) is subject to inventory constraints. Please inquire before ordering.

### ordering note

An extra charge is applied for cutting and/or straightening stainless steel and/or copper tubing, calculated from the total number of pieces produced for each line item

### cleaned

All tubing is cleaned and ready to use.

### Cleaned Copper Tubing

- Adheres to ASTM B-280.
- Chromatographically free of background contamination.
- Use for plumbing GC systems.

ID	OD	Wall	qty.	cat.#
0.065"	1/8"	0.030"	50 ft.	21590
0.190"	3/4"	0.030"	50 ft.	21592

### GC Installation Kit

This kit contains the tubing and fittings needed to add an additional GC to your lab bench. Also included are four 1/8-inch tees, so carrier, fuel, and other GC gases can be routed to the new inlet or detector from existing gas lines. Order additional parts, such as purifiers or regulators, separately to customize the GC installation to your specifications. Kit includes: one tubing cutter, one 1/8-inch x 1/4-inch reamer, one 7/16-inch wrench, one 1/2-inch wrench, four 1/8-inch brass tees, ten 1/8-inch brass nuts, ten brass front and back ferrules, and 50 feet (15.2 meters) of our instrument-grade cleaned 1/8-inch copper tubing.

Description	qty.	cat.#
GC Installation Kit	kit	21325

For the *ultimate*  
in inertness, use

**SILTEK®**  
tubing.

See pages 372-374.

### tech tip

#### Plumbing a GC

It is essential to use clean chromatographic-grade tubing to plumb a GC. Standard-grade tubing contains residual hydrocarbon contaminants from the drawing process. These contaminants migrate into the gas stream, elevating background noise and causing down time.

**1/16-Inch Tubing Cutter**

- Produces square, smooth cuts in 1/16-inch tubing.
- Eliminates tubing distortion.
- Replaceable cutting wheel.

Description	qty.	cat.#
1/16" Tubing Cutter	ea.	20192
Replacement Cutting Wheels	3-pk.	20185



**Ridgid™ Tubing Cutter**

- Excellent for cutting 1/8- or 1/4-inch metal tubing.
- Compact size is ideal for tight spaces.
- Replaceable cutting wheel.

Description	qty.	cat.#
Ridgid™ Tubing Cutter for 1/8" or 1/4" metal tubing	ea.	23011
Replacement Cutting Wheels	2-pk.	23012



**Tubing Reamer**

- Removes burrs from stainless steel tubing.
- For 1/4- or 1/8-inch tubing.
- Non-slip safety design.

Description	qty.	cat.#
Tubing Reamer	ea.	20134



**SSI TC-20 Tube Cutting Machine**

- Cuts 1/16", 1/8", or 1/4" tubing with inside diameter as small as 0.008".
- Electrically operated bench-top model.
- Handy dressing tool on the swing arm removes burrs and reams tubing.
- Voltage selectable 110–120/220–240 volts, 50–60Hz.\*

Description	qty.	cat.#
SSI Tubing Cutter Machine	ea.	23029
SSI Replacement Cutting Wheels	3-pk.	23030

\*Unit shipped set for 110–120 operating voltage. Switch to 220–240 volts by using alternate fuse and power cord (included).



CE

8" x 6 1/4" x 4 1/4"  
(20.3 x 15.9 x 10.8 cm)  
Weight: 11 lbs. (5.0 kg)

**Tubing Dressing Tool**

Same tool as included with the SSI tube cutting machine.

Description	qty.	cat.#
1/16" Tubing Dressing Tool	ea.	20188
Replacement Insert for 1/16" Tubing Dressing Tool	ea.	20189
1/8" Tubing Dressing Tool	ea.	20190
Replacement Insert for 1/8" Tubing Dressing Tool	ea.	20191



**Tubing Bender**

- Bends 1/8-inch, 3/16-inch, or 1/4-inch tubing.
- Accurate left-hand, right-hand, or offset bends.

Description	qty.	cat.#
Tubing Bender	ea.	23009



**1/16-Inch Tubing Cutting Pliers**

- Ideal for cutting 1/16-inch tubing.
- Cuts quickly, reducing distortion.
- Cuts clean, eliminating need for deburring.

Description	qty.	cat.#
1/16" Tubing Cutting Pliers	ea.	20193








## Tube & Pipe Fittings

### also available

For Silcosteel®, Siltek®/Sulfinert®, and Silcosteel®-CR-treated fittings, see pages 370-371.

### Tube & Pipe Fittings

The Instrumentation Group of Parker Corporation designs and manufactures a top-quality line of components and systems for use in process instrumentation, semiconductor manufacturing, and analytical equipment. Parker's product quality and delivery have made them a world-wide leader—and this is the level of quality and service Restek wants to deliver to you.

Fitting Type	Parker #	Similar to Swagelok® #	Size	qty.	Brass cat.#	316 Grade Stainless Steel	
						qty.	cat.#
nut 	1 Nu 1	102-1	1/16"	20-pk.	21800	5-pk.	21900
	2 Nu 2	202-1	1/8"	40-pk.	21801	10-pk.	21901
	4 Nu 4	402-1	1/4"	40-pk.	21802	10-pk.	21902
front ferrule 	1 FF 1	103-1	1/16"	20-pk.	21803	10-pk.	21903
	2 FF 2	203-1	1/8"	40-pk.	21804	20-pk.	21904
	4 FF 4	403-1	1/4"	40-pk.	21805	20-pk.	21905
back ferrule 	1 BF 1	104-1	1/16"	20-pk.	21806	10-pk.	21906
	2 BF 2	204-1	1/8"	40-pk.	21807	20-pk.	21907
	4 BF 4	404-1	1/4"	40-pk.	21808	20-pk.	21908
nut & ferrule set 	—	—	1/16"	10-pk.	21809	2-pk.	21909
	—	—	1/8"	20-pk.	21810	5-pk.	21910
	—	—	1/4"	20-pk.	21811	5-pk.	21911
plug 	1 BLP 1	100-P	1/16"	5-pk.	21815	2-pk.	21915
	2 BLP 2	200-P	1/8"	10-pk.	21816	4-pk.	21916
	4 BLP 4	400-P	1/4"	10-pk.	21817	4-pk.	21917
union 	1 SC 1	100-6	1/16"	3-pk.	21818	ea.	21918
	2 SC 2	200-6	1/8"	5-pk.	21819	2-pk.	21919
	4 SC 4	400-6	1/4"	5-pk.	21820	2-pk.	21920
reducing union 	2 RU 1	200-6-1	1/8" to 1/16"	5-pk.	21823	ea.	21923
	4 RU 1	400-6-1	1/4" to 1/16"	5-pk.	21824	2-pk.	21924
	4 RU 2	400-6-2	1/4" to 1/8"	5-pk.	21825	2-pk.	21925

Parker's (A-Lok) two-piece ferrules and NPT fittings are ideal for installing new equipment, modifying existing instrumentation, or replacing worn connections.

Size	Parker #	Similar to Swagelok®	qty.	Brass cat.#	316 Grade Stainless Steel		Fitting Type
					qty.	cat.#	
1/16"	1 ET 1	100-3	2-pk.	21826	ea.	21926	 tee
1/8"	2 ET 2	200-3	2-pk.	21827	ea.	21927	
1/4"	4 ET 4	400-3	2-pk.	21828	ea.	21928	
1/8"	2 ECR 2	200-4	2-pk.	21829	ea.	21929	 cross
1/4"	4 ECR 4	400-4	2-pk.	21830	ea.	21930	
A 1/8" tube to B 1/16"	2 TUR 1	100-R-2	5-pk.	21831	2-pk.	21931	 tube end reducer
1/4" tube to 1/16"	4 TUR 1	100-R-4	5-pk.	21832	2-pk.	21932	
1/8" tube to 1/4"	2 TUR 4	400-R-2	5-pk.	21833	2-pk.	21933	
1/4" tube to 1/8"	4 TUR 2	200-R-4	5-pk.	21834	2-pk.	21934	
1/8"	2 PC 2	201-PC	5-pk.	21835	2-pk.	21935	 port connector
1/4"	4 PC 4	401-PC	10-pk.	21836	2-pk.	21936	
1/8" tube to 1/4"	2 PC 4	401-PC-2	5-pk.	21837	2-pk.	21937	
A 1/8" tube to B 1/8" NPT	2 MSC 2N	200-1-2	10-pk.	21841	2-pk.	21941	 male connector
1/4" to 1/4" NPT	4 MSC 4N	400-1-4	10-pk.	21842	2-pk.	21942	
1/16" to 1/8" NPT	1 MSC 2N	100-1-2	5-pk.	21843	2-pk.	21943	
1/8" to 1/4" NPT	2 MSC 4N	200-1-4	10-pk.	21844	2-pk.	21944	
1/4" to 1/8" NPT	4 MSC 2N	400-1-2	10-pk.	21845	2-pk.	21945	
A 1/8" tube to B 1/8" NPT	2 FSC 2N	200-7-2	5-pk.	21846	2-pk.	21946	 female connector
1/4" to 1/4" NPT	4 FSC 4N	400-7-4	5-pk.	21847	2-pk.	21947	
1/4" to 1/8" NPT	4 FSC 2N	400-7-2	5-pk.	21848	2-pk.	21948	
1/8" male*	2A-Q4VN	QC4D-200	—	—	ea.	21957	 male & female quick couplings
1/8" male	2A-Q4P	QC4S-200	—	—	ea.	21958	
1/8" female*	2A-Q4CN	QC4B-200	—	—	ea.	21959	
1/4" male*	4A-Q4VN	QC4D-400	—	—	ea.	21960	
1/4" male	4A-Q4P	QC4S-400	—	—	ea.	21961	
1/4" female*	4A-Q4CN	QC4B-400	—	—	ea.	21962	

\*Includes self-sealing shut-off valve.

## Valves, Leak Detector, GC Installation Kit



Toggle-valve—instant on/off control



Ball valve—leak-free bidirectional sealing



Plug valve—leak-free at wide temperature swings



Metering valve with optional knob—accurate control of low flow



CE

Easy-to-clean probe assembly.



### Parker Shut-Off Valves

Parker toggle valves are ideal for applications in which instant on/off control is necessary. They are rated to 200psig at 21°C and have a maximum operating temperature of 148°C. Ball valves have a floating ball to assist sealing and to reduce operating torque, and dual seats to provide leak-tight bidirectional sealing. They are rated to 1500psig at 21°C and have a maximum operating temperature of 177°C. Perfect for instrument supply lines, plug valves work well in any application requiring throttling or on/off operations. Parker plug valves remain leak-free even when subjected to wide temperature swings. They are rated to 3000psig at 21°C and have a maximum operating temperature of 205°C.

Valve Type	1/8" Brass		1/4" Brass		1/2" Stainless Steel		3/4" Stainless Steel	
	cat.#	ea.	cat.#	ea.	cat.#	ea.	cat.#	ea.
Toggle	22188		22189		22190		22191	
Ball	22192		22193		22194		22195	
Plug	22196		22197		22198		22199	

### Parker Precision Metering Valves

Precision metering valves provide accurate, stable control of low gas and liquid flow rates. The valve stem threads do not contact the fluid stream, making these valves ideal for high-purity applications. The O-ring seal and stem threads are coated with a low vapor pressure, silicone lubricant for optimum performance. An optional vernier turn-counter knob allows repeatable flow settings for standardized operating processes.

Valve Type	1/8"		1/4"		1/2" Stainless Steel		3/4" Stainless Steel	
	Nickel-Plated Brass cat.#	ea.	Nickel-Plated Brass cat.#	ea.	cat.#	ea.	cat.#	ea.
Metering, straight	22200		22201		22204		22205	
Metering, right angle	22202		22203		22206		22207	

Vernier Knob for Metering Valve: cat.# 22209 (ea.)

### Restek Electronic Leak Detector

- Reliable thermal conductivity leak detector.
- Easy to hold and operate.
- Sensitive—detects helium at  $1.0 \times 10^{-5}$  cc/sec. or nitrogen at  $1.4 \times 10^{-3}$  cc/sec.
- Responds to leaks in less than 2 seconds.
- Audible alarm plus LED readout.
- Auto zeros with the touch of a button.
- Built-in rechargeable battery.

Description	qty.	cat.#
Leak Detector with 110Volt Battery Charger	ea.	22451
Leak Detector with 220Volt European Battery Charger	ea.	22451-EUR
Leak Detector with 220Volt UK Battery Charger	ea.	22451-UK

Caution: The Restek Electronic Leak Detector is NOT designed for determining leaks of combustible gases. A combustible gas detector should be used for determining combustible gas leaks under any condition. The Restek Electronic Leak Detector may be used for determining trace amounts of hydrogen in a GC environment only.

### Leak Detector Accessory Kit

The kit includes an adaptor fitting to detect leaks in hard-to-reach locations, and a mounting bracket that can be affixed to the wall or GC.

Description	qty.	cat.#
Leak Detector Accessory Kit (adaptor fitting for probe, mounting bracket)	kit	22453

### GC Installation Kit

This kit contains the tubing and fittings needed to add an additional GC to your lab bench. Also included are four 1/8-inch tees, so carrier, fuel, and other GC gases can be routed to the new inlet or detector from existing gas lines. Order additional parts, such as purifiers or regulators, separately to customize the GC installation to your specifications. Kit includes: one tubing cutter, one 1/8-inch x 1/4-inch reamer, one 7/16-inch wrench, one 1/2-inch wrench, four 1/8-inch brass tees, ten 1/8-inch brass nuts, ten brass front and back ferrules, and 50 feet (15.2 meters) of our instrument-grade cleaned 1/8-inch copper tubing.

Description	qty.	cat.#
GC Installation Kit	kit	21325

**54-Piece Tool Kit**

A must-have for every lab—all the tools you need in one place! Set comes with screwdrivers, pliers, wrenches, sockets, scissors, clamps, and more. Durable, zippered, notebook-style carrying case for easy transport.



Description	qty.	cat.#
54-Piece Tool Kit	kit	23004

**Plier Set**

Includes 6-inch nose/side cutter, 6-inch wire cutter, and 6-inch adjusting pliers.

Description	qty.	cat.#
Plier Set	set	23033

**Metric 9 Piece Ball-Point Hex Key Set**

Includes 9 metric hex keys (Allen wrenches): 1.5, 2, 2.5, 3, 4, 5, 6, 8, 10mm.

Description	qty.	cat.#
Metric 9 Piece Ball-Point Hex Key Set	set	22999

**12 Piece Ball-Point Hex Key Set**

Includes 12 hex keys (Allen wrenches): .050", 1/16", 5/16", 3/32", 7/64", 1/8", 9/64", 5/32", 3/16", 7/32", 1/4", and 5/16".

Description	qty.	cat.#
12 Piece Ball-Point Hex Key Set	set	22998

**Tubing Bender**

- Bends 1/8-inch, 3/16-inch, or 1/4-inch tubing.
- Accurate left-hand, right-hand, or offset bends.

Description	qty.	cat.#
Tubing Bender	ea.	23009

**Torx® Screwdriver Set**

- Set includes TR-10, TR-15, and TR-20.
- Ideal for performing routine maintenance on Agilent 6890 GCs.

Description	qty.	cat.#
Torx® Screwdriver Set	set	23034

**5-in-1 Magnetic Screwdriver**

Magnetic power tip holds bits and screws securely.

Description	qty.	cat.#
5-in-1 Magnetic Screwdriver	set	23002

**Ratchet Wrenches**

Easier to use in confined spaces, compared to adjustable wrenches.

Description	qty.	cat.#
3/8" Ratchet Wrench	ea.	23005
1/2" Ratchet Wrench	ea.	23006
7/16" Ratchet Wrench	ea.	23007
9/16" Ratchet Wrench	ea.	23008

**Wrench Set**

Includes 4-inch, 6-inch, and 8-inch adjustable wrenches.

Description	qty.	cat.#
Wrench Set	set	23001

