



Capillary GC Columns



Fused Silica MXT (Silcosteel)

**Widest Selection
of Columns
Available from ANY
Manufacturer**

***OVER 1200 Different
Sizes/Types***

IDs from 0.10, 0.18, 0.25, 0.32, 0.53mm

Length from 10, 20meter
15, 30, 60, 105meter

Film Thickness from
0.10, 0.20, 0.25, 0.50, 1.00, 1.50,
3.00, 5.00 (and 7.00um)

But depends on phase and ID
and column type
(see detailed catalog)

How do you sort ALL these possible types out ?

How do you optimise the column selection ?

What about maximising the performance of the selected column ?

***USE EzGC / Pro EzGC
Method Development Software***

**Restek
Australian
Distributors**

Chromalytic Technology Pty Ltd
Fax : +61 3 9761 1169 Phone : +61 3 9762 2034
E-mail : jimjeffs@chromtech.net.au

**Copyright Restek Corp.
and Chromtech : 1997
- All rights reserved**



Table of Contents for GC Columns by Phase

Custom Column Special Testing	18	Rtx-225	34
Integrated Guard Columns I 0	19	Stabilwax / MXT-WAX	35)
Capillary 6-Packs	19	Rtx-Wax	36
Rtx-1 / MXT-1	20	Stabilwax-DA	37
Rtx-1 PONA	20	Stabilwax-DB	38
Rtx-2887/MXT-2887	22	Rtx-2330	39
Rtx-5 / MXT-5	23	Rtx -502.2 / MXT-502.2	40
Rtx®-5MS	24	Rtx-624 / MXT-624	41
XTI-5	25	Rtx-Volatiles / MXT-Volatiles	42
Rtx-5 Amine	26	Rtx-BAC1 & BAC2	45
Rtx-20 / MXT-20	27	MXT-BAC1 & BAC2	43
Rtx-1301 / MXT-1301	28	Rt-bDEXm	44
Rtx-1701 / MXT-1701	29	Rt-TCEP	44
Rtx-35	30	Rt-Alumina	45
Rtx®-50	31	Rt-QPLOT / MXT-QPLOT Rt-S PLOT	45
Rix-200	32	Rt-Msieve	46
Rtx-65 / Rtx-65TG	33		

That sets Restek apart from the competition ?

Quality and Innovation you can count on!

High Quality Fused Silica Tubing

A columns quality is only as good as the materials from which it is manufactured. Restek uses only high quality fused silica tubing for our Rtx® capillary columns. We constantly monitor the strength of our tubing by using a manual winding system which assures that every meter of every column is checked. Our tubing suppliers often say, "If our tubing is strong enough for Restek, it's strong enough for anybody."

Innovative Cage Design

To keep your column safe we have designed a stainless steel cage to completely surround the column. The column is suspended in the cage using high temperature string that acts as a shock absorber. At no time does the fused silica tubing come into contact with the metal support of the cage. Spontaneous breakage can occur when fused silica tubing contacts sharp metal edges on the column cage. Our fused silica column cage is 7.65" (19.4 cm) in diameter and fits into virtually any GC oven.

Unbreakable Metal Capillary Columns

However, we realize that there are some harsh environments or special applications where fused silica tubing may not be the ideal material. For these circumstances Restek offers MXT capillary columns. MXT® columns are made from unbreakable, thin wall, stainless steel that has been treated with our Silcosteel process. This process deposits a submicron layer of flexible fused silica on the inside of the tubing giving it the same inertness as an Rtx® column, but with the ruggedness of stainless steel. MXT® columns are caged in small diameter coils that are ideal for portable or process control GCs.

Restek capillary columns are manufactured under US patent 4,293,415, licensed by Hewlett-Packard Company.



Restek
Australian
Distributors

Chromalytic Technology Pty Ltd
Fax : +61 3 9761 1169 Phone : +61 3 9762 2034
E-mail : jimjeffs@chromtech.net.au

Copyright Restek Corp.
and Chromtech : 1997
- All rights reserved



High Quality Polymers


In addition to the tubing material, high quality polymers used for deactivation and stationary phases is also necessary for making the world's best capillary columns. Most of the polymers that go into Restek's Rtx and MXT columns are synthesized to exacting standards in our own laboratory to ensure the highest quality and consistency. Any residual catalyst, which can cause degradation and increased bleed, is removed from the polymer. The polymer is then carefully fractionated to remove low molecular weight fragments. This provides a tight mono-modal distribution and further reduces column bleed. Every polymer is fully characterized to ensure long-term reproducibility. Characterization testing includes Refractive Index (RI), Fourier Transform Infrared Spectroscopy (FTIR), Kovats Retention Indices, % crosslinking, efficiency, and a five-day thermal bake-out to ensure column longevity. Each polymer lot is tracked in an extensive database for future reference and guarantees that the column you buy today will meet the same specifications of the column you bought last week.

Thorough Testing of Every Column

The final quality assurance test on every Restek capillary column confirms that each column has the necessary inertness and efficiency. Each column is also evaluated for bleed at its maximum operating temperature to ensure that every column exhibits the lowest bleed possible.

Functional, Environmentally Safe Packaging

The Restek difference is apparent before the column box is even opened. Instead of a logo or an advertisement on the top of the compact, recyclable corrugated cardboard column box, we've printed a useful Column Service Record. This allows the analyst to easily track column installation and GC maintenance for troubleshooting purposes. By having the column history on hand for instant review, quicker decisions can be made and guesswork eliminated.



*Neil Mosesman
Director of New Products &
Method Development*



*Rick Crago
Fused Silica
Production Supervisor*



*Kristi Sellers
FS Manufacturing Chemist*

Reduce, Resuse Recycle

In response to Restek's efforts to reduce and recycle, we were the proud recipients of the **1995 Governors Waste Minimization Award** sponsored by the Pennsylvania Department of Environmental Protection. This award recognizes companies for their accomplishments in preventing or reducing the waste in Pennsylvania.

Restek
Australian
Distributors

Chromalytic Technology Pty Ltd
Fax : +61 3 9761 1169 Phone : +61 3 9762 2034
E-mail : jimjeffs@chromtech.net.au

Copyright Restek Corp.
and **Chromtech** : 1997
- All rights reserved