

NitroGenerator



TEXOL
Products

Gas Generators

HROMalytic www.chromtech.net.au sales@chromtech.net.au

ABN 14 643 445 058 **PTY LTD** Tel: (03) 9762 2034 Fax: +61 3 9761 1169

Australian Distributors: **TECH**nology . . . for ALL Your chromatography supplies !

NitroGenerator

The award-winning NitroGenerator is a flexible and expandable state-of-the-art Nitrogen Generator that produces flows from 1 to 10 l/min at very high purity.




Overview

Winner of the "Millennium Products" award from the design council of Great Britain, the NitroGenerator is a discrete, efficient, modular system that produces high purity nitrogen only when you need it – all you need is compressed air. The NitroGenerator uses PSA technology for superior performance and reliability. A range of models is available to cover most applications. When selecting a system, consideration should be given to the required flow rate, purity and delivery pressure. The larger units are fitted with castors to allow them to be manoeuvred easily into position below a bench top. All of the units in the range are built with a tank module and control cabinet to allow the unit to be expanded upon.

Process Description

NitroGenerator has a built-in standby system, which ensures the most efficient use of electricity and the systems internal components. Compressed air enters the generator via a 0.01-micron filter [1] that removes contaminants and bulk moisture via an integrated auto drain. When the control system assesses that more nitrogen is required, solenoid valve [2] is opened and the pressure is allowed to build in the column [4] containing Carbon Molecular Sieve (CMS) and a residue moisture catcher [3]. This sieve adsorbs the oxygen and leaves the dry nitrogen free to flow into the reservoir [7]. After this point, the Nitrogen is allowed to flow to the application via pressure and flow regulators [8].



Features and Benefits

- Flow Rates from 1 to 10 litres/min
- High Purity nitrogen up to 99.9995%
- Very compact
- Light weight construction
- Stand-by Mode
- PSA technology
- Output Pressure 80 PSIG (5.4 Barg)
- Low power requirement
- Easy Installation
- Simple maintenance

Typical Applications

- Gas Chromatography (GC)
- Total Organic Carbon (TOC)
- Vessel Leak Detection
- Flame Ionisation Detector (FID)
- GC-MS
- Laser Beam Optics Purging
- ICP Spectrometer

Technical Specifications

Model:	NG1	NG2	NG3	NG5	NG10
Flow l/min	1	2	3	5	10
Purity (% N ₂)	Up to 99.9995% (5 ppm residual O ₂)				
Dewpoint (C PDP)	-40				
Inlet Pressure (Barg / PSig)	Purity dependent (7 Barg minimum)				
Outlet Pressure (Barg / PSig)	5.4 / 80				
Method	Carbon Molecular Sieve – Pressure Swing Adsorption				
Power (W)	8.1				
Voltage	Plug-in transformer (supplied) converts voltages from 110-230VAC to 24VDC				
Depth (mm / ins)	143 / 5.6	286 / 11.3	397 / 15.6	397 / 15.6	397 / 15.6
Height (mm / ins)	645 / 25.4	645 / 25.4	763 / 30.0	763 / 30.0	763 / 30.0
Width (mm / ins)	288 / 11.3	288 / 11.3	397 / 15.6	397 / 15.6	794 / 31.3
Weight (kg / lbs)	14.5 / 32.0	29 / 63.9	70 / 154.4	70 / 154.4	140 / 308.7
Connection Ports	Female ¼" Rc (PT) suitable for BSPT, PT, R & Unithread)				
Control Panel (h x w x d, wt)	280mm(11.4 ins) x 300mm(12.2 ins) x 125mm(5.1 ins), 3kg(6.6 lbs)				


HRMalytic www.chromtech.net.au sales@chromtech.net.au
 ABN 14 643 445 058 **PTY LTD** Tel: (03) 9762 2034 Fac: +613 9761 1169
 Australian Distributors: **TECH**nology . . . for ALL Your chromatography supplies !