


■ Syringes

■ Economical 700 Series (Hamilton)

- Cemented needles.
- Plungers and barrels are not interchangeable or replaceable.



700 Series - Solid Needle	Needle		Hamilton Cross Reference #	Restek Catalog #	Price
	gauge	style			
Hamilton 701 10ul LC Syringe for Rheodyne	22s	3	80365	21250	\$23.00
Hamilton 702 25ul LC Syringe for Rheodyne	22s	3	80465	21251	\$24.00
Hamilton 705 50ul LC Syringe for Rheodyne	22s	3	80565	21252	\$26.00
Hamilton 710 100ul LC Syringe for Rheodyne	22s	3	80665	21253	\$28.00
Hamilton 725 250ul LC Syringe for Rheodyne	22s	3	80765	212.54	\$34.00

■ Replacement Needles for Hamilton Syringes

Rheodyne style

22# gauge, 3 point style

for 10 to 100 µl syringe		
Hamilton Cross Reference #	Restek Cat. #	Price 3-pack
80464	21266	\$13.00
for 250 to 500 µl syringe		
Hamilton Cross Reference #	Restek Cat. #	Price 3-pack
80763	21267	\$13.00

Waters style

25s gauge, 3 point style

for 10 to 100 µl syringe		
Hamilton Cross Reference #	Restek Cat. #	Price 3-pack
80426	21268	\$13.00
for 250 to 500 µl syringe		
Hamilton Cross Reference #	Restek Cat. #	Price 3-pack
80726	21269	\$13.00


■ 1700 Gas Tight Series (Hamilton)

- Teflon Tipped plungers.
- Removable needles.
- Replaceable syringe barrels, plungers, and plunger tips.

1700 Series - Gas Tight Removable Needle	Needle		Hamilton Cross Reference #	Restek Catalog #	Price
	gauge	style			
Hamilton 1701 10ul LC Syringe for Rheodyne	22s	3	80065	21260	\$40.00
Hamilton 1702 25ul LC Syringe for Rheodyne	22s	3	80265	21261	\$42.00
Hamilton 1705 50ul LC Syringe for Rheodyne	22s	3	80965	21262	\$44.00
Hamilton 1710 100ul LC Syringe for Rheodyne	22s	3	81065	21263	\$46.00
Hamilton 1725 250ul LC Syringe for Rheodyne	22s	3	81165	21264	\$50.00

■ 800 Series (Hamilton)

- Reinforced plungers.
- Removable needles.
- Plunger/barrel assemblies are replaceable.
- Plunger stop prevents plunger blowout.



800 Series - Removable Needle	Needle		Hamilton Cross Reference #	Restek Catalog #	Price
	gauge	style			
Hamilton 801 10ul LC Syringe for Waters U6K	25s	3	8481.5	212.55	\$64.00
Hamilton 802 2ul LC Syringe for Waters U6K	25s	3	84816	21256	\$64.00
Hamilton 805 50ul LC Syringe for Waters U6K	25s	3	84817	21257	\$68.00
Hamilton 810 100ul LC Syringe for Waters U6K	25s	3	84818	21258	\$70.00
Hamilton 825 50ul LC Syringe for Waters U6K	25s	3	84819	21259	\$75.00

Prices are US List Prices in USD - Check with your local distributor for current prices

Restek : Australian Distributor - Chromalytic Technology Pty Ltd - Fax +61 3 9761 1169.