




Integrated Guard Columns/Column Six Packs

□ Innovative Integra-Guard Columns

How many times has it been said about press-fit connectors, "Some people swear by them. Others swear at them"? For many analysts, the aim of attaching a guard column to an analytical column is a mystery. Restek's chemists have discovered the solution to this mystery - the most reliable connection is no connection at all! No guard column system is more permanent than one continuous length of tubing containing both the guard column and the analytical column.



NEW !

*Guard columns
WITHOUT presstight
connections*

*protecting your
analytical column has
never been this easy !*

Restek now offers a wide variety of Integra-Guard™ capillary columns with a guaranteed leak-free connection! The column is tied using high temperature string with the guard column set apart from the analytical column. The transition area between the columns is the point at which the guard column ends and the analytical column begins. The entire setup is suspended in our unique "crush-free" cage that prevents the column from coming in contact with anything that could damage it.

The Integra-Guard™ column is so economical that we challenge you to beat our price with that of a conventional connection even if you assemble it yourself. If you are currently using a guard column or considering one for the future call Restek today and ask about our Integra-Guard™ columns.

mm ID	Length	Suffix #
0.25	5m	-124
	10m	-127
0.32	5m	-125
	10m	-128
0.53	5m	-126
	10m	-129

Ordering is easy. Just add the appropriate suffix number and price to the analytical column's catalog number and price.

For example. a 30m 0.25mm ID 0.25um Rtx®-5 with a 5-meter Integra-Guard column is cat.# 10223-124.



Jeanne Ladner
MXT Accessories
Group Leader

□ Spectacular Savings with Our Six Packs of Columns

Restek's most popular capillary columns are available in economical six-packs! Get six columns for the price of five. Other phases and configurations are available upon request.

Column Description	cat.#
30m, 0.25mm ID, 0.25um Rtx®-5	10223-600
30m, 0.25mm ID, 0.50um Rtx®-5	10238-600
30m, 0.32mm ID, 1.00um Rtx®-5	10254-600
30m, 0.53mm ID, 1.50um Rtx®-5	10270-600
30m, 0.32mm ID, 1.0um Stabilwax®	10654-600
30m, 0.53mm ID, 1.0um Stabilwax®	10655-600

You can
**SAVE
MONEY-**
get six columns
for the
price of five

Restek
Australian
Distributors

Chromalytic Technology Pty Ltd
Fax : +61 3 9761 1169 Phone : +61 3 9762 2034
E-mail : jimjeffs@chromtech.net.au

Copyright Restek Corp.
and **Chromtech** : 1997
- All rights reserved