

SRI Model 8610C

NEW for 1996



What's the Difference

New Chassis
Bright LED Display of All Operating Zones - states, temperatures
Digital Flow Control on ALL gases

The same Options and Accessories are Available as on the previous 8610 and 9300 models
PLUS !

Optional

Dual Oven - Get 2 GCs for the price of One Turbo-cool and sub-ambient temperature control.
Dual Oven Option
AND that's about half the cost of other GCs as well.

Main Features

Fast Cooling Oven - ambient to 400C
Electronic Pressure Control of all gases.
At-a-glance LED Display of ALL zones

Detectors

Choice of up to 4 detectors on each GC - FID, TCD, PID, ECD, NPD, FPD, and the Unique to SRI . . . Dry ELCD.
everything except MS (yet to come !)

Injection Systems

From 6 Different injection systems
1/8, 1/4" Packed Column (heated); Direct On-column (Cold or Heated); Split/Splitless
PLUS !
EPA Style Purge & Trap
Soil Desorption System
Gas Sample Valves (Heated and non-heated)
10 Port P&T Autosampler
AutoOnline Liquid Sampler
Headspace Accessory

**ALL COMPLETELY
CONTROLLED BY
THE PEAK SIMPLE
DATA SYSTEM**

Using either DOS or Windows